


# HPLC Method Description for In-process Control of Intermediates


<b>Project</b>	<b>Pyronaridine_INV-054926</b>
<b>Compound</b>	<b>PNDa01, PNDa02, PNDa06, PNDa04 and PND</b>
<b>Purpose</b>	<b>Method Description</b>
<b>Category</b>	<b>Methods</b>
<b>Substance Type</b>	<b>Intermediate</b>
<b>Report ID</b>	<b>INV_054926_HPLC_M1 Version 1.0</b>

## Authors

Cai Alain

## Approvers

By:  9DCE20AEFA87456 Title: AD Manager Date: 1/21/2024

By:  55362BFE31B84BD... Title: Head of Technical Development Date: 1/21/2024

## Distribution

Zhou Irvin, Pyronaridine (INV-054926) project team, Approvers

## Table of contents

1.	Objective	3
2.	Summary and conclusion	3
3.	Experimental	6
3.1.	HPLC	6
3.2.	Equipment and reagents	7
3.3.	Solutions	7
3.3.1.	Diluent	7
3.3.2.	Mobile phase preparation	7
3.3.3.	Solution preparations	7
3.4.	Proposed injection sequence and system suitability test	9
3.5.	Calculation	10
4.	Figures	11









## 1. Objective










This method 'INV\_054926\_HPLC\_M1' for intermediates of Pyronaridine (INV-054926) project is developed by HPLC. The parameters of the reversed phase HPLC method suitable for in-process control of intermediate steps shall be described in this document.




## 2. Summary and conclusion

The method is applied for the HPLC testing of all intermediates (in-process control). This method is based on reversed phase liquid chromatography with UV detection and gradient elution using a Waters Atlantis T3, 3µm, 150 x 4.6mm HPLC column.

**Table 1** Structure, Retention time and RRT of all intermediates and their related substances

Compound	Structure	Retention time (RT)	RRT
<b>PNDa01 step</b>			
6-methoxy-3-aminepyridine (SM 2)		ca. 6.6 min	ca. 0.34
PNDa01 impurity 2		ca. 17.4 min	ca. 0.91
2,4-dichlorobenzoic (SM 1)		ca. 18.0 min	ca. 0.94
PNDa01 impurity 1		ca. 18.2 min	ca. 0.95
PNDa01		ca. 19.2 min	1.00
PNDa01 impurity 3		ca. 20.5 min	ca. 1.07
<b>PNDa02 step</b>			
PNDa01		ca. 19.2 min	ca. 0.85
BIA		ca. 21.9 min	ca. 0.97

Compound	Structure	Retention time (RT)	RRT
PNDa02		ca. 22.5 min	1.00
<b>PNDa06-HCl step</b>			
PNDa06		ca. 7.7 min	1.00
Acetaminophen (SM5)		ca. 8.1 min	ca. 1.05
<b>PNDa04-HCl step</b>			
PNDa04		ca. 6.4 min	1.00
PNDa06		ca. 7.7 min	ca. 1.20
<b>PND step</b>			
PNDa04		ca. 6.4 min	ca. 0.54
DIA		ca. 11.4 min	ca. 0.97
PND		ca. 11.8 min	1.00
DIN		ca. 13.4 min	ca. 1.14

Compound	Structure	Retention time (RT)	RRT
TGF-001 impurity 1		ca. 16.2 min	1.37
TGF-001 impurity 2		ca. 17.1 min	1.45
PNDa02		ca. 22.5 min	1.91

Example chromatograms and extracted HPLC-PDA spectra are given in Section 4.

### 3. Experimental

Equivalent equipment or grade of materials can be used.

#### 3.1. HPLC

- HPLC System: Quaternary pump module (e.g.: Waters Alliance 2695)  
PDA detector (e.g.: Waters Alliance 2998)  
Auto sampler (e.g.: Waters Alliance 2695)  
Column oven (e.g.: Waters Alliance 2695)
- Empower-control and integration software or equivalent
- Column: Waters Atlantis T3, 3 $\mu$ m, 150 x 4.6mm
- Flow rate: 1.0 mL/min
- Elution: Gradient mode
- Run time: 30.0 min
- Detection:

PNDa01 step	254 nm
PNDa02 step	254 nm
PNDa06-HCl step	254 nm
PNDa04-HCl step	278 nm
PND step	278 nm

- Injection: 10  $\mu$ L
- Column temp.: 35°C  $\pm$  5°C
- Auto sampler temp.: Room temperature
- Mobile phase (see section 3.3.2):
  - **A:** 0.1% TFA in Water
  - **B:** 0.1% TFA in Acetonitrile
- Diluent:

PNDa01 step	1M Na <sub>2</sub> CO <sub>3</sub> in water/Acetonitrile=90/10 v/v
PNDa02 step	2-MeTHF
PNDa06-HCl step	Acetonitrile/ water = 20/80 v/v
PNDa04-HCl step	Acetonitrile/ water = 5/95 v/v
PND step	Acetonitrile/ water = 50/50 v/v

- Needle wash: Water/ Acetonitrile (50:50 v/v)
- Equilibration time: 5 min
- Gradient:

**Table 2** Gradient Table

Time (min)	% A	% B
0.0	100.0	0
3.0	100.0	0
13.0	65.0	35.0
20.0	1.0	99.0
25.0	1.0	99.0
25.1	100.0	0

Time (min)	% A	% B
30.0	100.0	0

### 3.2. Equipment and reagents

- Balance: e.g.: Mettler Toledo XP56
- Acetonitrile: HPLC grade, e.g.: Merck LiChrosolv
- Water: HPLC grade, e.g.: from Millipore ultra-pure water system
- TFA: HPLC grade, e.g.: Sigmer-Aldrich
- 2-MeTHF: HPLC grade, e.g.: Innochem
- DMSO: HPLC grade, e.g.: Innochem
- Glassware: 10, 50, 100 mL volumetric flasks, 1L graduated cylinders
- Pipette: e.g.: 1.0, 5.0 mL Pipette

### 3.3. Solutions

#### 3.3.1. Diluent

*Different volumes can be prepared as soon as the solvent ratio is the same.*

Diluent 1: 1M Na<sub>2</sub>CO<sub>3</sub> in water/Acetonitrile=90/10 v/v for PNDa01

Diluent 2: 2-MeTHF for PNDa02

Diluent 3: Acetonitrile/ water = 20/80 v/v for PNDa06-HCl

Diluent 4: Acetonitrile/ water = 5/95 v/v for PNDa04-HCl

Diluent 5: Acetonitrile/ water = 50/50 v/v for PND

#### 3.3.2. Mobile phase preparation

*Preparation is described for a volume of 1 liter. Different volumes can be prepared as soon as the solvent ratio is the same.*

Mobile phase A (0.1% TFA in Water):

In a suitable container, add 1000 mL of water and 1 mL of TFA. Mix well.

Mobile phase B (0.1% TFA in Acetonitrile):

In a suitable container, add 1000 mL of acetonitrile and 1 mL of TFA. Mix well.

#### 3.3.3. Solution preparations

*Other volumes and weigh-ins might be used as long as the final concentration remains the same. Min. weight of used balance must be considered during sample preparation.*

##### 3.3.3.1. Standard solutions

- PNDa01 step

IPC Standard Solution (10% SM1, 10% SM2 in 100% PNDa01):

Weigh approx. 2 mg of SM1, 2 mg of SM2 and 20 mg of PNDa01 reference standards into a 100-mL volumetric flask. Dissolve and dilute to volume with diluent 1. Mix well (SM1/ SM2: 0.02 mg/mL, PNDa01: 0.2 mg/mL).

- PNDa02 step

IPC Standard Solution (10% PNDa01 in 100% PNDa02):

Weigh approx. 0.8 mg of PNDa01 and 8 mg of PNDa02 reference standards into a 100-mL volumetric flask. Dissolve and dilute to volume with diluent 2. Mix well (PNDa01: 0.008 mg/mL, PNDa02: 0.08 mg/mL).

- PNDa06-HCl step

IPC Standard Solution (10% SM5 in 100% PNDa06):

Weigh approx. 3 mg of SM5 and 39 mg of PNDa06-HCl reference standards into a 100-mL volumetric flask. Dissolve and dilute to volume with diluent 3. Mix well (SM5: 0.03 mg/mL, PNDa06 free base: 0.3 mg/mL).

- PNDa04-HCl step

IPC Standard Solution (10% PNDa06 in 100% PNDa04):

Weigh approx. 6.5 mg of PNDa06-HCl and 83 mg of PNDa04-HCl reference standards into a 100-mL volumetric flask. Dissolve and dilute to volume with diluent 4. Mix well (PNDa06 free base: 0.05 mg/mL, PNDa04 free base: 0.5 mg/mL).

- PND step

IPC Standard Solution (10% PNDa02, 10% PNDa04 in 100% PND):

Weigh approx. 1.8 mg of PNDa02, 3 mg of PNDa04-HCl and 18 mg of PND reference standards into a 100-mL volumetric flask. Dissolve and dilute to volume with diluent 5. Mix well (PNDa04 free base/ PNDa02: 0.018 mg/mL, PND: 0.18 mg/mL).

### 3.3.3.2. Sample solutions

The IPC sample is submitted as a solution of each intermediate. Dilute the sample to a suitable concentration based on the estimated reaction solution concentration.

- PNDa01 step

For example: If the reaction solution is 360 mg/mL of PNDa01, transfer 0.5 mL of sample into a 50-mL volumetric flask and dilute to volume with diluent 1 (3.6 mg/mL), and then transfer 3 mL of that solution into a 50-mL volumetric flask and diluent to volume with diluent 1 (0.216 mg/mL).

- PNDa02 step

For example: If the reaction solution is 110 mg/mL of PNDa02, transfer 0.5 mL of sample into a 50-mL volumetric flask and dilute to volume with diluent 2 (1.1 mg/mL), and then transfer 0.7 mL of that solution into a 10-mL volumetric flask and diluent to volume with diluent 2 (0.077 mg/mL).

- PNDa06-HCl step

For example: If the reaction solution is 410 mg/mL of PNDa06 free base, transfer 0.5 mL of sample into a 50-mL volumetric flask and dilute to volume with diluent 3 (4.1 mg/mL), and then transfer 0.8 mL of that solution into a 10-mL volumetric flask and diluent to volume with diluent 3 (0.328 mg/mL of free base).

- PNDa04-HCl step



For example: If the reaction solution is 100 mg/mL of PNDa04 free base, transfer 0.5 mL of sample into a 50-mL volumetric flask and dilute to volume with diluent 4 (1 mg/mL), and then transfer 5 mL of that solution into a 10-mL volumetric flask and dilute to volume with diluent 4 (0.5 mg/mL of free base).

- PND step

For example: If the reaction solution is 230 mg/mL of PND, transfer 0.5 mL of sample into a 50-mL volumetric flask and dilute to volume with diluent 5 (2.3 mg/mL), and then transfer 0.8 mL of that solution into a 10-mL volumetric flask and dilute to volume with diluent 5 (0.184 mg/mL).

### 3.4. Proposed injection sequence and system suitability test

**Table 3** Proposed injection sequence and SST criteria

Sample name	No. of injections	SST acceptance criteria
Blank (sample diluent)	1 + N <sup>[1]</sup>	No interference between the blank peaks and the components of interest
IPC standard Solution	1	For identity
Sample solution	1	N/A

<sup>[1]</sup> Additional blanks may be run until an acceptable baseline is obtained.

### 3.5. Calculation

Reaction monitoring by formula (%area):

- **PNDa01 step**

$$\% \text{ SM2 (remaining by conversion)} = \frac{A_{SM2}}{A_{SM2} + A_{PNDa01}} * 100\%$$

Where:

$A_{SM2}$  = peak area of SM2 in sample solution (diluted reaction solution)

$A_{PNDa01}$  = peak area of PNDa01 in sample solution (diluted reaction solution)

- **PNDa02 step**

$$\% \text{ PNDa01 (remaining by conversion)} = \frac{A_{PNDa01}}{A_{PNDa01} + A_{PNDa02}} * 100\%$$

Where:

$A_{PNDa01}$  = peak area of PNDa01 in sample solution (diluted reaction solution)

$A_{PNDa02}$  = peak area of PNDa02 in sample solution (diluted reaction solution)

- **PNDa06-HCl step**

$$\% \text{ SM5 (remaining by conversion)} = \frac{A_{SM5}}{A_{SM5} + A_{PNDa06}} * 100\%$$

Where:

$A_{SM5}$  = peak area of SM5 in sample solution (diluted reaction solution)

$A_{PNDa06}$  = peak area of PNDa06 in sample solution (diluted reaction solution)

- **PNDa04-HCl step**

$$\% \text{ PNDa06 (remaining by conversion)} = \frac{A_{PNDa06}}{A_{PNDa06} + A_{PNDa04}} * 100\%$$

Where:

$A_{PNDa06}$  = peak area of PNDa06 in sample solution (diluted reaction solution)

$A_{PNDa04}$  = peak area of PNDa04 in sample solution (diluted reaction solution)

- **PND step**

$$\% \text{ PNDa02 (remaining by conversion)} = \frac{A_{PNDa02}}{A_{PNDa02} + A_{PND}} * 100\%$$


Where:

$A_{PNDa02}$  = peak area of PNDa02 in sample solution (diluted reaction solution)


$A_{PND}$  = peak area of PND in sample solution (diluted reaction solution)

#### 4. Figures


**Figure 1** Example HPLC chromatogram of PNDa01 IPC solution, method INV\_054926\_HPLC\_M1, 254nm




**Figure 2** Example HPLC chromatogram of PNDa02 IPC solution, method INV\_054926\_HPLC\_M1, 254nm



**Figure 3** Example HPLC chromatogram of PNDa06-HCl IPC solution method INV\_054926\_HPLC\_M1, 254nm



**Figure 4** Example HPLC chromatogram of PNDa04-HCl IPC solution, method INV\_054926\_HPLC\_M1, 278nm



**Figure 5** Example HPLC chromatogram of PND IPC solution, method INV\_054926\_HPLC\_M1, 278nm

